

Edition 10/2013
Replace Ed. 04/2008



Soy Peptone A2 SC 19649

Organotechnie® S.A.S.

27, avenue Jean Mermoz
93120 La Courneuve, France
Tél : +33 (0) 1 49 92 87 50
Fax : +33 (0) 1 49 92 87 51

e-mail : info@organotechnie.com
web : www.organotechnie.com

Definition

Soy Peptone A2 SC is manufactured by a controlled enzymatic hydrolysis of soybean meal.

Description

Fine beige powder easily soluble in water.

Soy Peptone A2 SC contains a mix of peptides and free amino acids.

Soy Peptone A2 SC is manufactured with raw materials from vegetal origin only.

Soy Peptone A2 SC is guaranteed "non GMO" according to the European Regulations 1829/2003 and 1830/2003

Use

Source of organic nitrogen recommended in media for:

- Analytical microbiology
- Industrial fermentation.

Physico-chemical characteristics

	Standard
Solubility in water at 5 %	Complete
pH (5 % solution)	6.4 - 7.8
Loss on drying	≤ 6 %
Total nitrogen TN	8.5 - 11.0 %
α-amino nitrogen AN	2.5 - 3.5 %
AN/TN x 100	23 - 41
Residue on ignition	≤ 20 %
Chloride (as NaCl)	≤ 1 %

Microbiology

	Standard
Total aerobic microbial count	≤ 10 000 /g
Coliforms	≤ 10 /g
<i>Escherichia coli</i>	Absence / g
<i>Salmonella</i>	Absence / 25 g
<i>Staphylococcus aureus</i>	Absence / 10 g
Yeasts and moulds	≤ 20 /g



The information contained in this publication is based on our own research and development work and is to the best of our knowledge true and accurate.

Users should, however, conduct their own tests to determine the suitability of our products for their own specific purposes.

Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for the infringement of any patents.

Soy Peptone A2 SC / 19649

Edition 10/2013
Replace Ed. 04/2008




Organotechnie S.A.S.

Typical data


Molecular weight distribution	g / 100 g
> 10 000 daltons	0
1 000 - 10 000 daltons	8.6
500 - 1 000 daltons	25.9
300 - 500 daltons	34.9
< 300 daltons	30.6

Average Molecular Weight **503 daltons**



Molecular weight distribution

Amino Acids		Total - T (g/100 g)	Free - F (g/100 g)	F/T x 100
Alanine	Ala	2.3	0.3	13.0
Arginine	Arg	4.3	1.6	37.2
Aspartic acid	Asp	6.0	0.1	1.7
Cystine	Cys	2.1	0.2	9.5
Glutamic acid	Glu	10.2	0.4	3.9
Glycine	Gly	2.4	0.5	20.8
Histidine	His	1.3	0.2	15.4
Isoleucine	Ileu	2.1	0.2	9.5
Leucine	Leu	3.9	1.1	28.2
Lysine	Lys	3.6	1.0	27.8
Methionine	Met	0.8	0.2	25.0
Phenylalanine	Phe	2.1	0.4	19.0
Proline	Pro	2.5	0.1	4.0
Serine	Ser	2.9	0.8	27.6
Threonine	Thr	2.2	0.2	9.1
Tryptophan	Trp	0.6	0.2	33.3
Tyrosine	Tyr	2.0	0.4	20.0
Valine	Val	2.5	0.3	12.0



Amino Acids F/T x 100

Documentation

The certificate of analysis and the sanitary certificate are supplied with each delivery.

Packing and storage

25 kg net corrugated board box with inner polyethylene bags.

Upon request: 5 kg plastic drum.

Keep in original packaging closed when not in use,
at room temperature in a dry area.

Hygroscopic product.

Best before: 5 years.

Health and safety information

Dusty powder.

Avoid inhalation.

The information contained in this publication is based on our own research and development work and is to the best of our knowledge true and accurate.

Users should, however, conduct their own tests to determine the suitability of our products for their own specific purposes.

Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for the infringement of any patents.